



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C11orf57-chromosome 11 open reading frame 57 Gene

Reference: PTXBC030546

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf57; chromosome 11 open reading frame 57		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 57
GENE ID	55216
ORIGIN SPECIES	Human
ACCESSIONS	BC030546
VECTOR	pDONR223
NCBI SYMBOL	C11orf57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagacagatggggfcacagtggttataaagagttataccctgaagaatttgaacacagacagatgagatcatcagaagatattaccaacgggaagaaaaatctccccagtgaaagtctatcccatgaatcccgaacacaaagaagtcaaagaaatccccacaaaaaagcagaaa
aaaaggtcacacaaaaaacagaagaaaagcaaaaaaggaagccacagatatatacagcagattctctgagtgagttctcagaanaactggggctctctgtacaaaggaaggaacacacataaacccaagaaaaaatccaggaaaaaagctctcataaaaactctgttattcttag
ggcagaagaatcaacttcacatcagatgattcagcatcagatcgtctctgagtgaggaagagacactaagaaaaccaaaggaaaaaggaaggaagaaaaagccatactctgtgccaacaatgaataacagagagaggacaaacacgcacaattggaaagtgtctac
acagatgaaaagctctctcgtcagtcagagatgactaa