

C10orf118-chromosome 10 open reading frame 118 Gene

Reference: PTXBC030557

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf118; coiled-coil domain-containing protein 186; CTCL tumor antigen HD-CL-01/L14-2; CTCL tumor antigen L14-2; coiled-coil domain containing 186		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 118
GENE ID	55088
ORIGIN SPECIES	Human
ACCESSIONS	BC030557
VECTOR	pDONR223
NCBI SYMBOL	C10orf118
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagagacagaccatagcctctactctctgataaaaatgttgggaaacacctaattaaagggaagactcatgcaactgttttctggcaatgaaagcagcaaatagaaaacgagtcctcaactattgtcattaaacactgataaaactttatgtcaacctaagagcataataatcgaa
ttgaagcccaggaaaattatattccagatcatgttgagggtgaggattctgtgcccacacagacaggtctcagaaaattctgaacaaatagctaatctctagtggaattttgctaaacatatattcaaaaacaaatgaacagacagaaaagtaacacaaatattggtggaattaaggatc
ctacatttccagaatcagctaatgaaaagacttattcagaagccctatgatacagactgcaccaagaaatttcaaaaaataaagagcgtttcagcatcagaggattgttggagaataagaatctgagctcttatctacggagtttgagaacatcgagtaccaaattggaatgaataagg
gagaacatgcattagtgtgttgaaggtgtgtgcaagataaatatttcagcaggaacatacataaaaaagtaa