

LOC376693-hypothetical LOC376693 Gene

Reference: PTXBC030568

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC376693
GENE ID	376693
ORIGIN SPECIES	Human
ACCESSIONS	BC030568
VECTOR	pDONR223
NCBI SYMBOL	LOC376693
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgatgcctaagaagaactggattgccatttatgaactccttttaaggaggaggatcatgggtggccaagaaggatgtccacatgcctaagcacccggagctggcagacaagaatgcgccaaccttcgtgcatgaaggccatgcagtctctcaagtcgccgaggtacgtgaagggaagt
ttgcctggagacatttctactggtaccttaccatga