

RPS15A-ribosomal protein S15a Gene

Reference: PTXBC030569

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S15a; 40S ribosomal protein S15a; up-regulated by HBV X protein; ribosomal protein S15a		

For research use only.

GENE DESCRIPTION	ribosomal protein S15a
GENE ID	6210
ORIGIN SPECIES	Human
ACCESSIONS	BC030569
VECTOR	pDONR223
NCBI SYMBOL	RPS15A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcgcataatgtcctggcagatgctctcaagatfatcaacaatgccgaaaagagaggcaaacgccaggtgcttattaggccgtgtccaaagtcacgtccggtttctcactgtgatgaagcatggttacattggcgaattgaaatcattgatgaccacagagctgggaaattg
ttgtgaacctcacaggcaggctaaacaagtgtgggtgatcagccccagatttgacgtgcaactcaaaagacctggaaaaatggcagaataatctgttccatcccgcagtttggtttcattgtactgacaacctcagctggcatcatggacatgaagaagcaagacgaaaacacacagga
gggaaatcctgggattcttttctag