

DKKL1-dickkopf-like 1 (soggy) Gene

Reference: PTXBC030581

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT34; SGY; SGY-1; SGY1; dickkopf-like protein 1; cancer/testis antigen 34; dickkopf-like 1; protein soggy-1; soggy; soggy-1; testicular secretory protein Li 15; dickkopf like acrosomal protein 1		

For research use only.

GENE DESCRIPTION	dickkopf-like 1 (soggy)
GENE ID	27120
ORIGIN SPECIES	Human
ACCESSIONS	BC030581
VECTOR	pDONR223
NCBI SYMBOL	DKKL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagaagcctcccacctgccccgcaaggcggcatctgctggtctctgctgctctctctacccctggtgatccccctcctgctgagctcctatccatgatctgacgccaagagagctcttgggtctcacaggcctccagagcctactccaaggctcagccgactttctgaaa
ggtaacctgctcggggcatagacagcttattctgccccatggacttcggggcctcctgggaactaccacaaagagagagaccaggagcaccagctggggaacaacacctctccagccacctccagatcgacaagatgaccgacaacaagacaggagaggtgctgatctccg
agaatgtgtggcatcattcaaccagcggaggggagctcaggggtgattgaagggtaccaggatggaggagaaggaggccctggtaccatccagaaggccacggacagcttccacacagaactccatccccgggtggcctctggatcattagctgccacggcggaggtccc
accaggatgccctggaggcggccactggctcagcgagaagcgacaccgctgcaggccatccgggatgactccgcaaggggaccacaaggacgtctagaagaggggacgagagctcctccactccaggctgtccccgaaagacccacttactgtatcctcaggcc
ctctcgagctgtag