

PPP3R2-protein phosphatase 3 (formerly 2B), regulatory subunit B, beta isoform Gene

Reference: PTXBC030595

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPP3RL; calcineurin subunit B type 2; CBLP; CNBII; calcineurin B, type II (19kDa); calcineurin B-like protein; calcineurin BII; protein phosphatase 2B regulatory subunit 2; protein phosphatase 3 (formerly 2B), regulatory subunit B (19kD), beta isoform (ca		

For research use only.

GENE DESCRIPTION	protein phosphatase 3 (formerly 2B), regulatory subunit B, beta isoform
GENE ID	5535
ORIGIN SPECIES	Human
ACCESSIONS	BC030595
VECTOR	pDONR223
NCBI SYMBOL	PPP3R2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccacaatgggaacgaggccagttacccggcggagatgtgtctccactttgacaatgatgaattaaaaggctgggcaggagggttaagaagttggacttgacaaatcagggtctctgagcgtggaggagttcatgtccctgccgagctgcgccacaacccgttggtgcggcga
gtgatcgacgtcttcgacaccgacgtgatggagaagtggaactcaaggaaattcatcctggggacctcccagttcagcgtcaaggcgcgagaggagcagaagttgaggtttgcgttcagcatttacgacatggataaagatggctacattccaacggggagctctccaggtctgaagat
gatgggtgggcaacaacctgacggactggcagctccagcagctggtcgacaaaacctcatcatcctggacaagatggcgtgggaagataatcctttgaggaaattcagtgctgtgtcagagacctggagatccacaagaagctgtgtctcatcgtatga