



Reference: PTXBC030603

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bA787I22.1; armadillo repeat-containing protein 2; armadillo repeat containing 2		
For research use only.			

GENE DESCRIPTION	armadillo repeat containing 2
GENE ID	84071
ORIGIN SPECIES	Human
ACCESSIONS	BC030603
VECTOR	pDONR223
NCBI SYMBOL	ARMC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr