



Reference: PTXBC030608

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCAN; SLRR5A; podocan proteoglycan		
For research use only.			

GENE DESCRIPTION	podocan
GENE ID	127435
ORIGIN SPECIES	Human
ACCESSIONS	BC030608
VECTOR	pDONR223
NCBI SYMBOL	PODN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr