

ATXN7L1-ataxin 7-like 1 Gene

Reference: PTXBC030616

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATXN7L4; ataxin-7-like protein 1; ataxin 7-like 4; ataxin-7-like protein 4; ataxin 7 like 1		

For research use only.

GENE DESCRIPTION	ataxin 7-like 1
GENE ID	222255
ORIGIN SPECIES	Human
ACCESSIONS	BC030616
VECTOR	pDONR223
NCBI SYMBOL	ATXN7L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgtcggagcgttctcgaatccgtgtctcggctgctgctgccgaaggaacagggaaaaagcaacaagaaggagcaatggcgacactggatcgcaaaagtgccagtcggaggcgtttctgggcaaacctggctcctggatcgacgccgcaaatcactgtcgcg
acaatgtagattagaagaggctgaaaaagagggtggaaaaagcaggaggttatgaggcttaataaagaagatagcacttattggccattaccagcacatgacgacttctctcgtagtgtgcagtgccctgtaaccaggtcgtaagccacagggtttccagtcgcactgcgggagaa
agcaagacaacaggagaaatgaaggcatctccaggagtgaccagagagcagccaagccatagagaagcatcagggtgtga