

RNF182-ring finger protein 182 Gene

Reference: PTXBC030666

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF182; ring finger protein 182		
For research use only.			

GENE DESCRIPTION	ring finger protein 182
GENE ID	221687
ORIGIN SPECIES	Human
ACCESSIONS	BC030666
VECTOR	pDONR223
NCBI SYMBOL	RNF182
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagtcacacctctgaagacactcggagctcagggcctctgatgagctggagtgcaaaatctgttacaatcgatacaatctgaacagaggaaacccaagtcctgagtggtgtcataggggtttgtgccaatgcctctacaagatcatagactttgggactccccacaaggtgc
attgtctgtcttctgcaggtttgagacgtgctgccagatgatgaagtagtagcctgccgatgacaacaacatcctgttaacttgacttggaggcaaaagggaagagtcctgccagagaacccctactgagctgctcaccaccaagaggtggcctctctgtcagtccttctca
cacgtcctccaaactgctggtcataaccatcatggaggtgcagagagagctccccgtccctgagctccactcctgtgtgagaattttataggcctgcgagtttcgactctgtcaccactgtgtcacacaactggactgtgtggaactgcacgtccctgctgttcagacatccatccgggtgtt
agtgtggttgctaggtttgctctactcagctccttacccttaggaatctactactggtgtctaagaaagtacaccttggggctgctttgtcagcctgggtcccttcgagcctcgttattcttatgtgtatggttttccagtggtgtgtcatgaatttcagactgtatggcacctcctttaa