

TMEM148-transmembrane protein 148 Gene

Reference: PTXBC030801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM148; NCRNA00311; long intergenic non-protein coding RNA 311		
For research use only.			

GENE DESCRIPTION	transmembrane protein 148
GENE ID	197196
ORIGIN SPECIES	Human
ACCESSIONS	BC030801
VECTOR	pDONR223
NCBI SYMBOL	TMEM148
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggctgctgaccagggtagaggctggaacccctggatggcccaacctggaggtagctagccatgcctctgccattggccggccacacaggttctgccctcttctcgcgtgccctgtgcctcccgctgctcgcgtcctaaccatgctggactggcagagactgaaaa
cgccccacggggccctgtgctgctcactcgggctcctccggagagtggtggcaacccctgcgcctgagcatcttgctctgcttggtgcagctggatcccttttggacttggggctgcacacgcacctgtgtctggagcgaatgggtcgtgagagccgagccctga