

TMEM156-transmembrane protein 156 Gene

Reference: PTXBC030803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 156		
For research use only.			

GENE DESCRIPTION	transmembrane protein 156
GENE ID	80008
ORIGIN SPECIES	Human
ACCESSIONS	BC030803
VECTOR	pDONR223
NCBI SYMBOL	TMEM156
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaaaacagccctcctaaattattgtgcaatagtgatcacattttaaatttgcgggaattttcaagacaccgaaagaaagacattggagctatcatgtctggaagtgtgttgaacttaattttacctattcactctctccttaattttctttgtgactttctgcaaccagtaagg
gaaactcagattatcatgagaatttttaaatccctcaattttcgttaactcaccaggacttgccaagacatcacaggtaatttaaatgtgtctctctgtgttgggttgagcctaaaggaaacatggattttatttctcaggaaacacatcaaaagtcttatcaggagagatcaatggaag
tgaagcaaatgattttcattcactgtcagcactttaactcagtgtagctctctgtgtgaccactggagggaatataacactacctgtcatctaaaaaacacactggaagatcaacatcatggagatgagccaagcaaggagaaatcgataaactacactgtagaatcatggaatacc
cgaatgattgtatcacatttttgcacctggagatggatataaaaaatcactgttccatgaagatcacttggtatatttttagttctattagttttatattttgatcatcctcactatccgcaaaatacttgaaggccagagaagagtgcaaaagtggcagagtcataagacaaacctacatctgt
tctcttaagagggaagtgttcggagaactgagagcattgaatgtgcagggtctttcagcagagaccacgcagaggtcgtctttggatcaagtcagggaagtgtcttcccccaattccagaactataa