

IGKV1-5-immunoglobulin kappa variable 1-5 Gene

Reference: PTXBC030814

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	immunoglobulin kappa variable 1-5
GENE ID	28299
ORIGIN SPECIES	Human
ACCESSIONS	BC030814
VECTOR	pDONR223
NCBI SYMBOL	IGKV1-5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctccctgctcagctcctgggctgctaagctctgggtctctggatcagtggggatalgtgatgaactcagctcctcctgcccgtcaccctggagagccggcctccatcctgcaggtctagtcagagcctcctgcatagtgatggatacaactattggattgttacctg
cagaagccaggggcagctctccacagctcctgatctattgggttctaatacgggctcgggggtccctgacaggttcagtgccagtgatcagggcacagattttacactgaaaatcagcaaatggagctgaggtggtgatttactgcatgcaaggctacaactcctcagacgttcg
gccaaggaccaaaggtggaatcaaacgaactgtggctgcacacatctgtcttcatcttcccgcctctgatgagcagttgaaatctggaactgcctctgtgtgtgcctgctgaataacttctatcccagagagggccaaagtacagtggaaggtggataacacctccaatcgggtaactccc
aggagagtggtcacagagcaggacagcaaggacagcactacagcctcagcagcaccctgacgtgagcgaagcagactacgagaacacaaagtctacgctgcgaagtcacccatcaggggcctgagctcggcgtcacaaagagcttcaacaggggagagtgtag