

HPCAL4-hippocalcin like 4 Gene

Reference: PTXBC030827

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	hippocalcin-like protein 4; hippocalcin like 4		
For research use only.			

GENE DESCRIPTION	hippocalcin like 4
GENE ID	51440
ORIGIN SPECIES	Human
ACCESSIONS	BC030827
VECTOR	pDONR223
NCBI SYMBOL	HPCAL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaagaccaagcaagctggcccccagggtgctggaggacattgtcagaacactgagttcagcagcaggagctgaagcagtggtacaaagggtctgaaggactgccccagcggcctcctcaacctggaggagtttcagcagctctacatcaagttctccctacggcga
cgctccaagttcgcgagcagcgtttccgacattcgacaagaacggcgacggcaccatcgactccgggagttcatctcgccctgtcgggtacctcccgcgagcttcgagcagaagctcaactggcctttgagatgtacgacctggacggcgacggcgcatcagcgccctg
gagatgctggagatcatcgaggcaatctacaagatggtggcaccgtgatcatgatgcgcatgaaccaggacgggctcacgccccagcagcgtgtggacaagatctcaagaagatggaccaggataaggacgaccagattacattggaggagttcaaggaggcagccaagatgac
gcatccattgtgtgctgctgacgtgtgacatgcagaagtag