

CD83-CD83 molecule Gene

Reference: PTXBC030830

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD83 molecule; CD83 antigen (activated B lymphocytes, immunoglobulin superfamily); CD83 antigen; BL11; HB15; B-cell activation protein; cell surface protein HB15; cell-surface glycoprotein; hCD83		

For research use only.

GENE DESCRIPTION	CD83 molecule
GENE ID	9308
ORIGIN SPECIES	Human
ACCESSIONS	BC030830
VECTOR	pDONR223
NCBI SYMBOL	CD83
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcgcggcctccagcttctgctcctgagctgcgcctacagcctggctcccgcgacccggaggtgaaggtgcttgcctcgaagatgtggaactgccctgcaccgccccctgggatccgcaggttcctacacggctctgggtcaagttattggagggtgtgaagagaggatg
gagacacccaggaagaccacctcaggggacagcactatcatcagaaggggcaaaatggttcttgcacgcccccaatgaaaggccctattccctgaagatccgaaacactaccagctgcaactcgggacatacagtgcaactctcaggaccggatgggacagagaacctaagt
ggcaagtgatctgagagtgacaggatgccctgcacagcgtaaagaagagactttaagaaatacagagcggagattgctcctgctgctgctggttatttctacttaacacatcatcttcaactgtaagtttgacggctacagagatctccagattttctaaagctggcatggaacg
agcttttcccagttacctccccaaataagcatttagggctagtgactcctcacaagacagaactggtatga