



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SUCNR1-succinate receptor 1 Gene

Reference: PTXBC030948

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

GPR91; succinate receptor 1; G-protein coupled receptor 91; P2Y purinoceptor 1

For research use only.

GENE DESCRIPTION	succinate receptor 1
GENE ID	56670
ORIGIN SPECIES	Human
ACCESSIONS	BC030948
VECTOR	pDONR223
NCBI SYMBOL	SUCNR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgggggagcatfagcgaatgcaacttgcgaaaactgctgctgcagcagagctgcccctgaaaaagtactactcttccatttttatgggattgagtcgttggggagctcttggaaataccattgtgtttacggctacatctctctctgaagaactggaacagcagtaatatattatctttta
accctctctctgactgactgtttcttgcacccctcccatgctataaggagttatgccaatggcaaaactggatataggagacgtctctgtaagaacacccgatatgtctgcatgcacacctctatcacgactctcttccacttttaccgcatagatgcatacttaataaagtatctct
ccggaacacctctctcaaaagaagaagtgtctatttaattctctggccacttgggtttttttatgaactctgactactaccccttataaalctcttttaactgacatgcaaaccccttgatattgttgcagctctggagaccccactacacactcatttaccagatgctgtctaac
actgttgggggtctctatctctttttgtgattgtttctttattacaagaatgctctctcttaaaagcagaggaaataggcaggttgcctactgctctgccccttgaaaagcctctcaactgtgctcatatggcagtggttaatactctctctgtttttacacccctatcagctatgcgggaatgtgagatgc
ctcacccggctggggatttgggaagcagatcatgctgactcagctctctatcaactctttatcatgtgcacacggccttggccctctctgaacagctgtcaacccctgctctcttatttcttttggagatcacttcaggagacgtctgatgaatcaactgagacacaaactcaaatcccttatcatcttt
ctcagaatggctctatgacactctctattcttcagaaaaagtga