

PRAC-prostate cancer susceptibility candidate Gene

Reference: PTXBC030950

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	small nuclear protein PRAC; PRAC; C17orf92; small nuclear protein PRAC1; prostate cancer susceptibility candidate protein 1; prostate, rectum and colon expressed gene protein; prostate cancer susceptibility candidate 1		

For research use only.

GENE DESCRIPTION	prostate cancer susceptibility candidate
GENE ID	84366
ORIGIN SPECIES	Human
ACCESSIONS	BC030950
VECTOR	pDONR223
NCBI SYMBOL	PRAC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgtgcgccatttctcagatcaaggaccggcccatcttactacctccaagagtgcttttctcttaataagaaaacatctactttgaacatctactggcgagaccaggagtgatggctcagcctgttaattctggaattcgggagggccgaggcaggaagattcctga