

TMEM135-transmembrane protein 135 Gene

Reference: PTXBC030952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PMP52; transmembrane protein 135; peroxisomal membrane protein 52		
For research use only.			

GENE DESCRIPTION	transmembrane protein 135
GENE ID	65084
ORIGIN SPECIES	Human
ACCESSIONS	BC030952
VECTOR	pDONR223
NCBI SYMBOL	TMEM135
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccctcagcaagtcacccatcataactgctatgagatcgccacacttggcaccctctgccgggtctctctcgcagatcaccggggggcgccctggaggagtcctgaagatctatgctcctctgtacttgattgcagcaattctccggaaacggaaattagactattattaca
caaaactacccctgagatcctacaatccgcttcatttctaaactgctaatggggccttgatattgctttttgcattttaaggagatacttgaaaaattctactcatggactcctggctttgggtgccgctctgccagcattctatgtggccattctcattgaaagaaaaagcaggagaggctgctc
acaatttatatggccaacttgatggagtctcactgtctcccccaggtggagtgatttcagctccctgcagcctctgcctccaggttcaaggagattctcctgcctga