



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ANK1-ankyrin 1, erythrocytic Gene

Reference: PTXBC030957

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANK; SPH1; SPH2; ankyrin-1; ANK-1; ankyrin 1, erythrocytic; ankyrin-R; erythrocyte ankyrin; ankyrin 1		

For research use only.

GENE DESCRIPTION	ankyrin 1, erythrocytic
GENE ID	286
ORIGIN SPECIES	Human
ACCESSIONS	BC030957
VECTOR	pDONR223
NCBI SYMBOL	ANK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgactcttcgtacccagctgtttgctacgcgtggtgctgtcagcttctcttgctcagctgcagaaagctgatgcacattctcaggggctcctgtgcttggctgaagacatccaccaggagctggacaaggagctggggagagcaggacctcagtgacgacgaggagacca
tctccacagggttgctgcggcgccgggtctcttcgaaggggaatgattgttcagaattatccaggggcagctgacagaggcaatccaggcggatgcacggcgcaacattgcaccagaagaatcattcgaagggtgtgtgcacagatagctgtccacgcgcgatgccgccaggca
gcacgcgaaggtgtgagctgcggagcgaatggtccctacacgcgcgcagctgacacggcgcaagaaagctgcagatagtagaagccgcgcacctcgaacaaagggtggaagaaacagtga