

SLC22A2-solute carrier family 22 (organic cation transporter), member 2 Gene

Reference: PTXBC030978

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OCT2; solute carrier family 22 member 2; organic cation transporter 2; solute carrier family 22 (organic cation transporter), member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 22 (organic cation transporter), member 2
GENE ID	6582
ORIGIN SPECIES	Human
ACCESSIONS	BC030978
VECTOR	pDONR223
NCBI SYMBOL	SLC22A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaccacgtggacgatgtcctggagcatggagggtttcacttttccagaagcaaatgttttctcttggctctgtctcggctaccttcgcgccatctacgtgggcatcgttctctgggttcaccctgaccaccgtgccggagccccggagtggccgagctgagctctgc
gctgcggctggagctcctgcagagggaactgaactacacgggtgccccggccagacgtgcgggagcctcccaagacagtgtaggcgctacgaggtggactggaaccagagcaccttcgactgcgtggacccccctggccagcctggacaccaagagccgctgccactgg
gccccctggcgacggcgtgggtgtacgagacgcctggctcgtccatcgtcaccgagtttaacctggatgtgccaactcctggatgttgacattccagtcacgtgaatgtaggattctttattggctctatgagatcggctacatagcagacaggtttggccgtaagctctgctcctaa
ctacagtcctcataaatgtcgacgtggagttctcatggccatttcccaacctatcgtggatgttaatttttcgcttaaccaagactgtcagcaaaagcaggctggtaaataggctacatcctgagtaagaatgtttgtcgtgcaactgtgagaacaaagcaacctcgtgccaaaatag