

## C11orf17-chromosome 11 open reading frame 17 Gene

Reference: PTXBC030996

|                |                                                                                                                                                                                                                                                     |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C11orf17; A-kinase-interacting protein 1; A kinase (PRKA) interacting protein 1; breast cancer associated gene 3; koyt binding protein 1; koyt binding protein 2; koyt binding protein 3; proline-rich protein BCA3; A-kinase interacting protein 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 11 open reading frame 17                                                               |
| <b>GENE ID</b>          | 56672                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC030996                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C11orf17                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggacaactgtttggcgccgcagcgcgtgaatgggtggaccgacgttccctgcagcgttcagcaaggctggcctagaagtgtggagaggcccaaggaggcgggtggactggcatgccctggagcgtcccaaggcgtgcattgggggtccttggccgggagggcgccacacct  
agagaaacagccggcagccggcccgagcgcgttctccgggagagagagaagagagaccccccaacccttagtgccttcagacaatggctgaattcatggactatactcaagtcagtggtggaaatattatctatctgtgccagagggaaggaggcgcaaccatgtctatcgttat  
cacagaggcgagtcgaagctgcacatgtgcttgacatagggaatggtcagagaaaaagacagaaaaagacatcccttggtcctggaggcagctatcaaatatcagagcatgctccagaggcatccagcctgctgagaacatcttaaggacctctacatagaagtatatccaggac  
ctattctgtcactgtggcctcaaatgacttaaccaagaagactcatgtggtagcagttgattctggacaaagcgtggacgtggtcttccctgtgtga