

MPHOSPH6-M-phase phosphoprotein 6 Gene

Reference: PTXBC031017

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPP; MPP-6; MPP6; M-phase phosphoprotein 6		
For research use only.			

GENE DESCRIPTION	M-phase phosphoprotein 6
GENE ID	10200
ORIGIN SPECIES	Human
ACCESSIONS	BC031017
VECTOR	pDONR223
NCBI SYMBOL	MPHOSPH6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgagagaaagacaaagtgtccaagaatctgtgcgcatgaagtgtatgcaaggaggactggactcagaacacagaacaaactagaagaagaagaaagaaatcattagtgaagagcactgggtactggattgccagagcttaagagaaagagagtttcataatagaa
gagcagagtttcttactatgtgaagatctctctatggaagaatgtcattcagaggatgtaacctgaggttgagaaatgatgcttcagatgaatgtaagcacaagcagaagaagtgaagatgaaacagtagagcttgatgtgtcagatgaagagatggctagaagatatgagaccttggtg
gggacaattgggaaaaagttgccagaagagagaccatgccaatatgaagaggatgaaatggagacataacacaaattaaagcaagaagatgttcttaagccccaggattaa