



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IFI30-interferon, gamma-inducible protein 30 Gene

Reference: PTXBC031020

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFI30, lysosomal thiol reductase; IFI-30; IP-30; IP30; gamma-interferon-inducible lysosomal thiol reductase; gamma-interferon-inducible protein IP-30; interferon gamma-inducible protein 30 preproprotein; interferon, gamma-inducible protein 30; legumaturai		

For research use only.

GENE DESCRIPTION	interferon, gamma-inducible protein 30
GENE ID	10437
ORIGIN SPECIES	Human
ACCESSIONS	BC031020
VECTOR	pDONR223
NCBI SYMBOL	IFI30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgtcgccactctctgctgttctgccaccgtctgtctgtctgtggacgtccccacggcggggtgcaggcgccccctctcaagcgttagactctttgggaatgggccaccagttaactacaagacaggcaatctatacctggggggccccctgaagaagtccaatgcaccgtgtgcaatgtgacctctactatgaagacgtgtgcgtggctgcgaacgtctctctgacccggagctcttcccaacatggctgttggtcatggagatctcctaatgtcacgtgtgtgccccacggaacacgacgaagcaaaaatgtcagtggcaggtgggaagtctcaatgtccagcatgggagaagaggagtgcaaaatcaacaaggtggagggctgcgtgttgatgaactgacatggagctagccctctgacactgtgtcgtatgaagagtlttagaacatggagagaagtctgccatagccctgcagctctacgccccaggcgctgcgcagacatacatgaggtgtgtcaaaagggggacccgggatgcagctatgcacgccacgccccggcgacagatgctctccagccaccgcagatgattgtgcctgtggcctgcctcaatggaaaccccttgaagatcagaccagctcttaccctgtgtccagtgtgtaccaggcggaagaccggatgtgtcccttctcaaccagctccctcaggaggtgtgcttccaagtga

More info: Please contact us for any question: contact@proteogenix.fr