

EPB41L4A-erythrocyte membrane protein band 4.1 like 4A Gene

Reference: PTXBC031042

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPB41L4; NBL4; band 4.1-like protein 4A; erythrocyte protein band 4.1-like 4; erythrocyte membrane protein band 4.1 like 4A		

For research use only.

GENE DESCRIPTION	erythrocyte membrane protein band 4.1 like 4A
GENE ID	64097
ORIGIN SPECIES	Human
ACCESSIONS	BC031042
VECTOR	pDONR223
NCBI SYMBOL	EPB41L4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctgtttctgcgtgtccggaagaattttactgcgaagtttgcctcctggatgaatccaagttaacccctaccaccagcagcagggcatcaagaagtcacgaaagggttcgttgccttgaccacgtattccatcacgtaaacctgtggagatagatttttgggctactgtactgtga
cagaagccatcagacgtattgctgacccctgcaaaaaccctgtgaacacaaagaactgatcaacactggaccccatatactttgtatttgggtattaaattctatgctgaagatccatgtaaacctaaagaagaataaccagatatacgtttttcttcaggtgaagcaagatgtccttcagg
gccgtctgccctgtcccgtcaacactgctgctcagctgggagcgtatgccatccagtcggagcttggagattatgacccatataaacatactgcaggatattgtctgagtagccggttgcctctgatcagaaggagaacttgaagaagccatagaaggattcataaaactctaattgtaa