

## PARP11-poly (ADP-ribose) polymerase family, member 11 Gene

Reference: PTXBC031073

|                |                                                                                                                                                         |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ARTD11; C12orf6; MIB006; poly [ADP-ribose] polymerase 11; ADP-ribosyltransferase diphtheria toxin-like 11; poly(ADP-ribose) polymerase family member 11 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | poly (ADP-ribose) polymerase family, member 11                                                    |
| <b>GENE ID</b>          | 57097                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC031073                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PARP11                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtttcacaagcagaagaattattttctaaacaacaacaatgaagtggatgacatggacacgtcagataccagtggtgctgttttacttggcagaatgtgggaagtggcagatgttcagccggataccaacagtcagtggtcagttagcagtgagatatacgaataaaagcttcaaaa  
caaacctgtgtgctcatttctttactacttcaaatcagctacaagatagactttgcagaaatgaagcaaatgaatctcaccactggaaagcagcgttaataaaaaagagccccctttctatcagtgctttcagttacatctgtgaaaacgaggccatccctatgccaccactgggagaat  
gtgaatactcaagtaacatcagcttattctctgcacaatcaaacatgaatataatgaagtgtccttcttgggaagacgatggatcgaaccgaataaaagaattcagagaattcaaacctagatttggagggtcttttgcaggaataaaggctcagctcaagaaaaaagaggt  
gtgcctcagattaatgaacaatgctgtttcagtggtaccagcagtgaaattgtggaagcaatctgcattcataactttgattggagaataaatgtatatacgtgtgtcttggaaaagataacatgtggaagctgtag