

PXT1-peroxisomal, testis specific 1 Gene

Reference: PTXBC031105

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	STEPP; peroxisomal testis-specific protein 1; small testis-specific peroxisomal protein; peroxisomal, testis specific 1		

For research use only.

GENE DESCRIPTION	peroxisomal, testis specific 1
GENE ID	222659
ORIGIN SPECIES	Human
ACCESSIONS	BC031105
VECTOR	pDONR223
NCBI SYMBOL	PXT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagctgagacattgggacacattgatcataggatgttcgagaggatctcaacaggatggcagagatgcactagatcatttgtctcttttctttaagaggtcaggtgttgcatttttctggaacaaccatttgcgttaa