



Reference: PTXBC031213

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ISL1 transcription factor, LIM/homeodomain; ISLET1; Isl-1; insulin gene enhancer protein ISL-1; islet-1; ISL LIM homeobox 1		

For research use only.

GENE DESCRIPTION	ISL LIM homeobox 1
GENE ID	3670
ORIGIN SPECIES	Human
ACCESSIONS	BC031213
VECTOR	pDONR223
NCBI SYMBOL	ISL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]