

FAM71C-family with sequence similarity 71, member C Gene

Reference: PTXBC031221

Product type: DNA Storage Condition: -20°C
Size/Quantity: 2 ug Delivery Condition: +4°C, Blue Ice
Synonyms: protein FAM71C; family with sequence similarity 71 member C
For research use only.

GENE DESCRIPTION	family with sequence similarity 71, member C
GENE ID	196472
ORIGIN SPECIES	Human
ACCESSIONS	BC031221
VECTOR	pDONR223
NCBI SYMBOL	FAM71C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggagactctgtatgctaccgtattacacggcccaagcagcccccaatggcatgttaacacctccatggggaagctgcagcggcaactgtacaaggcagtagatactatcaggtatgcaccatgtttgagagcactttatcagatcagcaaaaggagagaagtattga
cgtgcacaaccgtcccgaatggtacaatgggcatcgttcgaccagccccctgcctcacactctgatgcatgctgctggcccaccagctgctgctgacaaatgccaggtgtgctcctgcccaccagaaaagagaagtccgctcagagatctagaactaacaggctgctt
ccctgatgtttgtgaaataaccatccacaacagcgtaaaaaacagctccacctgaagcttgcacggcgcctctttatctcagctgtgtcctccatcggatgcaagtgaaagactttgttccactgggaaaacctgtttacatcctgagaccaccagtgaggcttacagtgatacca
ggcctatcctagctgggaacacattggactcatctgtttggaggaagtgcagaggagcccagtgggatgcccagaaagttttggaagagaaggagcaatcagaatcagtagactcacaatgctgaaatgtttgggtccacctattgtattacaatagatatga
```