



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

C17orf74-chromosome 17 open reading frame 74 Gene

Reference: PTXBC031286

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C17orf74; chromosome 17 open reading frame 74		
For research use only.			

GENE DESCRIPTION	chromosome 17 open reading frame 74
GENE ID	201243
ORIGIN SPECIES	Human
ACCESSIONS	BC031286
VECTOR	pDONR223
NCBI SYMBOL	C17orf74
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr