

ANKRD39-ankyrin repeat domain 39 Gene

Reference: PTXBC031303

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 39; ankyrin repeat domain 39		
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 39
GENE ID	51239
ORIGIN SPECIES	Human
ACCESSIONS	BC031303
VECTOR	pDONR223
NCBI SYMBOL	ANKRD39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgagcctcggccctgcgcggagggccctgctgctgcatccagcgggtgctcgggtacagcagacgctggaggagatggactcggaggggaatctgctcggcagccctgaatggagacctgggccgagtgaagcatttaacagaagggcggagaccaagtc
agcccgactcggcggtacactgcctgcatatgccagccgaatgggcactacgctgtgctccagttcctgctggaagcgggagctaatgtgatgccagaccacgggggtgacctgctgctgaccgagccagctactgcgggacactgaaatcacgcggctcctgctatc
acatgggtccaacccaggggtggtggtgacgacggcatgaccagttctgataaggctgctgagaggggtcacgggacatctgctccctcctcctgcaacacagccagccctgaaggccatccgggaccgaaaggcacgggtagcatgtgacctgtgcttgaacagtgacct
gcgggacctgtatccagctga