

LY6D-lymphocyte antigen 6 complex, locus D Gene

Reference: PTXBC031330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Ly-6D; lymphocyte antigen 6D; e48 antigen; lymphocyte antigen 6 complex, locus D		
For research use only.			

GENE DESCRIPTION	lymphocyte antigen 6 complex, locus D
GENE ID	8581
ORIGIN SPECIES	Human
ACCESSIONS	BC031330
VECTOR	pDONR223
NCBI SYMBOL	LY6D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggacagcattgtctccttgcaacctggctgtggtctacagggccagcccttacctgcgtgccagtggtgaccagctccagcaactgcaagcattctgtgtctgccggccagctctcgtctgcaagaccacgaacagtgaggcctctgagggggaatctggtgaag
aaggactgtgcggagtggtgcacaccagctacacctgcaaggccaggtcagcagcggcaccagctccaccagtgtctgcccagaggacgtgtgcaatgagaagctgcacaacgtgcacccaccgcaccgcctcgcaccagtgccctcagcctggggctggccctgagcc
tcttgccgctccttagccccagcctgtga