

POLE4-polymerase (DNA-directed), epsilon 4 (p12 subunit) Gene

Reference: PTXBC031331

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YHHQ1; p12; DNA polymerase epsilon subunit 4; DNA polymerase II subunit 4; DNA polymerase epsilon p12 subunit; DNA polymerase epsilon subunit p12; polymerase (DNA) epsilon 4, accessory subunit; polymerase (DNA-directed), epsilon 4 (p12 subunit); polymerase		

For research use only.

GENE DESCRIPTION	polymerase (DNA-directed), epsilon 4 (p12 subunit)
GENE ID	56655
ORIGIN SPECIES	Human
ACCESSIONS	BC031331
VECTOR	pDONR223
NCBI SYMBOL	POLE4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgggcgggcgggcagggaagcgggacgccccgagaggaggaggtacctgctggggaggcagcggcctcgagccccaggccccaacgagtgctgctgggctcgtctcaggttgccctctggcgcgagtgaggccttggtgaaggcagatcccgcgtgacgct
agcgggacaggaagccatcttctggtcacgagcgcggaactgtttgtggagaccattgcaaaagatgcctactgttgcgctcagcagggaagaaagcaaaccttcagaggagagacttgataatgcaatagaagctgtggatgaatttctttctggaaggtacttttagattga