

MRPS15-mitochondrial ribosomal protein S15 Gene

Reference: PTXBC031336

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DC37; MPR-S15; RPMS15; S15mt; 28S ribosomal protein S15, mitochondrial; MRP-S15; mitochondrial ribosomal protein S15		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S15
GENE ID	64960
ORIGIN SPECIES	Human
ACCESSIONS	BC031336
VECTOR	pDONR223
NCBI SYMBOL	MRPS15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgaggggtcgcgtggaggacgtgagtttgattcgaccggcgagttaccaggtcctagtagccgggtgcggcggtggagcgccaagtcttcaaccagtggggcctgcagcctcgaagtctcctccaggccgcgcggatgtcgtccggaaccagccc
agtcaggctggatgacccacctctctacgtgctcaaaactaccagaatgtccctggaattgagaagggttgatgagtcgtgaaaagactctgtcttggaaatggccaacaagaaggagatgctaaaaatcaagcaagaacagttatgaagaagattgttgcacccagagg
acaccagatccctggaggctcgaattattgcttctgtcaagatccgcagttatgaagaacacttggagaacatcgaaaggacaaagcccaaacgctatctgctaatgagcattgaccagaggaagaaagatgctcaaaaacctccgtaacaccaatgatgtcttgagaagatat
gctgggggtgggaattgagtacacctccccctctgtattaccgaagagccaccgaggtcgtgaccaagaaggctctgtgattcgggtttccaggagactcaaaagctgaagaagcgaagaagagccttaaggctgcagcagcagcccaaaacaagcaagcggaggaa
cccagacagccctgccaagccataccaagacactcaagacagccaataa