

C2orf58-chromosome 2 open reading frame 58 Gene

Reference: PTXBC031410

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf58; CYP1B1 antisense RNA 1		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 58
GENE ID	285154
ORIGIN SPECIES	Human
ACCESSIONS	BC031410
VECTOR	pDONR223
NCBI SYMBOL	C2orf58
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggctgttaataagctgagcttgaggcttcaaggactctgtgagctgtgctggatctagcagtgaccagtgatgccagtgccatggatgccagggccatgttcagtagtgagcaggtgctagcagcagaagcaatagaatggaaactgggagatctgctgtggatgaatcctca
ttatctctttctctgacctctttggactaaataaactccagttattcttgaccaacagaggtagggttggacagaacctcaaaagtcattagacttctactaataaggcattggaagtcaagcgactaactaggaaacaggttagttattgtaccctacactgtgcaacagattacatac
cattgtctatttggcctccttcttaggctagattccagagaagaagattactga