

ZFAND3-zinc finger, AN1-type domain 3 Gene

Reference: PTXBC031481

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TEX27; AN1-type zinc finger protein 3; testis expressed sequence 27; zinc finger, AN1-type domain 3; zinc finger AN1-type containing 3		

For research use only.

GENE DESCRIPTION	zinc finger, AN1-type domain 3
GENE ID	60685
ORIGIN SPECIES	Human
ACCESSIONS	BC031481
VECTOR	pDONR223
NCBI SYMBOL	ZFAND3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagacgtgggagcgcagcagcaaaagcgccagcctgcccgcctcgtccctgcggtctctgggggtccagcaagactatgaatctctgttccaaatgcttgcgtatttcaaaagaacagcagcagatgattccgctccaagtacaagtaacagccaatcagattgtttcc
gaagagaccaccagtgacaacaatactcgtataaccagccaactcttagtcccagccagcagccgttccgacagaactgaatgtaactccagcagtaaaagaggagtggtggccatgcacagacacagctcatgtctcattatcacaccaaaaaagatcctgtgtacagatt
cacagctctgagaatgaggcttcacagtaaaacggccagactacttgagaatacgggaacgggtccgaggaacacagtcgatctaaacagaagagtcgacgtcgggtgcttccagtgccaaacaaactggagctggtgcagcagggaattgggatcgtgtcgtcgtggttatgtttctgtat
gttacatgcctccccgagcagcagctgcacattcgaccacatgggcccgtggccgggaggaagccatcatgaaatggtgaagctggaccgaaagtggggcgctcctgccagcgcacgcgggaggggtgtcctga