

TRA@-T cell receptor alpha locus Gene

Reference: PTXBC031551

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRA@; IMD7; TCRA; TCRD; TRAC; T-cell antigen receptor, alpha polypeptide; T-cell receptor, alpha (V,D,J,C); T-cell receptor alpha locus		

For research use only.

GENE DESCRIPTION	T cell receptor alpha locus
GENE ID	6955
ORIGIN SPECIES	Human
ACCESSIONS	BC031551
VECTOR	pDONR223
NCBI SYMBOL	TRA@
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcctgaaattctcgtgtccattcttggattcagttggcatgggtgagcaccagctgtggagcagagccctcagtttctaaagcatcaagagggagaaaatctcactgtgtactgcaactcctcaagtgtttttccagcttacaatggtacagacaggagcctggggaagtcctgtc
ctcctgtgtgacagtagttacgggtggagaagtgaagaagctgaagagactaaccttctcagtttgggtgatgcaagaaaggacagttctctccacatcactgcggccagcctgtgtgatacaggcctctgtgtctaaacttggaaatgagaaattaacctttgggactggaacaagactcaccatc
atacccaatccagaaccctgaccctgcccgtgtaccagctgagagactctaaatccagtgacaagctgtctgtcctattcaccgattttgattctcaacaataatgtgtcacaaagtaaggattctgtgtatcagacacacacaaactgtgtacatgaggtctatggacttcaagagcaac
agtgtgtggcctggagcaacaatctgacttgcattgtcacaacgccttcaacaacagcattatccagaagacacctcttccccagccagaaaagttcctgtgtatgtcaagctgtgtcgagaaaagctttgaaacagatacgaactaaacttcaaacctgtcagtgattgggtccgaat
cctctcctgaaagtggcgggtttaatctgtctatgacgctgcggctgtgtgtccagctga