

C11orf75-chromosome 11 open reading frame 75 Gene

Reference: PTXBC031564

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0443 protein C11orf75; C11orf75; FN5; single-pass membrane and coiled-coil domain-containing protein 4; single-pass membrane protein with coiled-coil domains 4		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 75
GENE ID	56935
ORIGIN SPECIES	Human
ACCESSIONS	BC031564
VECTOR	pDONR223
NCBI SYMBOL	C11orf75
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccagctcaaggggaagcccaagaaggagacctccaaggacaagaaggagcgggaagcaagccatgcaggagcccgccagcagatcactacagtgggtactgccacgctggccgtggtcgtgctctgatcgtggtgtttgtgtacgtggccacgcgccccaccatcacga
gtga