

**FAM19A4-family with sequence similarity 19 (chemokine (C-C motif)-like), member A4
Gene**

Reference: PTXBC031566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM19A4; TAFA-4; TAFA4; chemokine-like protein TAFA-4; family with sequence similarity 19 (chemokine (C-C motif)-like), member A4; family with sequence similarity 19 member A4, C-C motif chemokine like		

For research use only.

GENE DESCRIPTION	family with sequence similarity 19 (chemokine (C-C motif)-like), member A4
GENE ID	151647
ORIGIN SPECIES	Human
ACCESSIONS	BC031566
VECTOR	pDONR223
NCBI SYMBOL	FAM19A4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagggtcccaaggatgagagctgtgctaaagtcaagtgctgctgcactggctcttctagcctacgtgttaatgggtgctgtaagctgatgctcgctcaagccagcacctccgggacatgaggtcaccaccaaatcaagcaaggacctgtgaggtggtcgcctgacacaggtgctgcaataagaaccgcatagaagagcgggtcacaacgggtcaagtgctcttctccgggacaggtggcgggcacaactcgggtcaacctctgtgtgaagcttcattgtgattcagaatgggtggtcaccatgaatccgtgtttggaaggagagattgaaagtgtccagattctcaggttgctctgtagcagtgaggcaataaagtcaaaactacgaaggtaacgcggtag