

MS4A15-membrane-spanning 4-domains, subfamily A, member 15 Gene

Reference: PTXBC031610

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	membrane-spanning 4-domains subfamily A member 15; membrane spanning 4-domains A15		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 15
GENE ID	219995
ORIGIN SPECIES	Human
ACCESSIONS	BC031610
VECTOR	pDONR223
NCBI SYMBOL	MS4A15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttcgccgcggccacgtggcgcctctctcatcgaggcggcgctccctctctggggaggagcctgcttcacatctcggatccctctcagtgccagccgagaagaaccacacagttgctggtgaggagcagcctgggcaccaacatcctcagcgtcatggcgccctttgctggg
acagccattctgctcatggatttgggtgtaccaacgggatgtggacagggcctatctggcgtgcttactatcttcaactgtcctggagttcttcacagcggtcattgccatgcacttgggtgccaagccatccatgccagggccagtgacacctgtgatcttctgccaacgccttcagcgca
gacttcaacatccccagccggcagcctctgcgccccctgcctatgacaatgtggcatatgcccaaggagtcgtctga