



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

KIAA1257-KIAA1257 Gene

Reference: PTXBC031632

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA1257		
For research use only.			

GENE DESCRIPTION	KIAA1257
GENE ID	57501
ORIGIN SPECIES	Human
ACCESSIONS	BC031632
VECTOR	pDONR223
NCBI SYMBOL	KIAA1257
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcgtctttttaccacattgagtattctcttcgccgacgatgaagaacctaaaaagttgacatattgtatttccaatgggtgccaaagtattctcggagtcaggagtaaaagactgtgaagcgttgccacgaagggtgacaagcctgggtgtcgtgggagcagacttttaatactactgtg
acaaaggaaatttataagaataaaatcccaaaaacaccttgaggctcttgaacactaaagaacagatgtcaagaaaagtcagatattaccgattaaagactcggcgcttcacagacgacgtgggagctttctataagtcaagaatgtagacactttgttttaatacagagaaaattatctgaa
caggggcatgtagaatacacaactgtcagacagagatgtcgaacaggaaacatccgccgaagaacaagaaaaacagaaaaacccaagattgttccaagttcttcacattacagagagcggccgaacttcttcagaagaacatgtaggaatatggaagatccctcaaaatggacgattcttc
cagcattcaatggaggtgtttcagaacacacacacattctttgcaggagcagatcagatggagactcaaggaaataattgagatgtaactcattgcgtttcaaacattacagacacaaaggtctcaaaatlaagggtggaagattgcaggaggaagaaatccagagagagacata
aagaagccctgtgcagaagaagagcagacaccagctgcagacgcccgaagcagagcgtctttcatcagcttgcgcgtcttcctgcgcgttcgaacactcgttgcctcttccagacactctcgtgtgtgtgtgtgcagacagcgcacataa

More info: Please contact us for any question: contact@proteogenix.fr