



Reference: PTXBC031653

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uracil nucleotide/cysteinyl leukotriene receptor; P2Y-like receptor; R12; UDP/CysLT receptor; G protein-coupled receptor 17		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 17
GENE ID	2840
ORIGIN SPECIES	Human
ACCESSIONS	BC031653
VECTOR	pDONR223
NCBI SYMBOL	GPR17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1