

C12orf54-chromosome 12 open reading frame 54 Gene

Reference: PTXBC031670

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C12orf54; HSD-29; HSD-30; chromosome 12 open reading frame 54		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 54
GENE ID	121273
ORIGIN SPECIES	Human
ACCESSIONS	BC031670
VECTOR	pDONR223
NCBI SYMBOL	C12orf54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacagcatccctgccaggatcaggaaacaaaggtagaatgacatccaagcagcagagagaagcacatccatagaagagacaatgagaccacaggaacacaggttaaccatcactgaaaccctgtggaccagggtgctgacagtttttaaggatatacaaaaggagctgcaggaa
gatgctcggattcgaggatgagcaactgctccatgacacccatgacatcagcaccaggactggaagcataaggcctccagattccttgatgaccccaagtggagaagattgcagttcagctctggagagcagccatcaggaggccgtatccacaacctgaagacacagctcttcagt
caatcagcttactaccctggaccctaa