

## TEKT3-tektin 3 Gene

Reference: PTXBC031688

|                               |                                                             |                     |                |
|-------------------------------|-------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                         | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | tektin-3; testicular microtubules-related protein; tektin 3 |                     |                |
| <b>For research use only.</b> |                                                             |                     |                |

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tektin 3                                                                                         |
| <b>GENE ID</b>          | 64518                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC031688                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TEKT3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggaacgtgtaggtgtactttaacgacaacttacgccacccctagaccaacacccaactttctaccagccatgaccatggcctcaagctacaggagcgtttcccccactccaatttgaccatagcctgagcctctctggagaccagcacatactacaagtgcctccaattc  
cccaagcgtggcccgctactgaccagatcacagaggggtgctcgagaataccatgcttccctttgttccacagaaacacactttctcaagatacacaccggatgactgggtacaggtccaatttaaccaactatcaagagtcacaacttcccgacataatcggagaaactaagatggat  
acatctcgctgattcaagacaaatatacaacaagaaaaactcaggcagacacacccaaatctgggagacgtgtcaatgacataggggttttggaaatctgaaatcattcatgagttggatgaaatgattggagagacaaatgcacttactgatgtgaagaaaagactggagcgggc  
tttgatggagactgaagccctctctcaggtagcccgagaatgtctatttctcagaaaaagagaatgggaatcgacatgctcagatgaagtgaagcacaactgctgacgggaagtgtgatactattctgtgtgtcaagaaagaatgaagctacatttggataaggctattgcccaacttgacg  
ccaacagagcgtccagcatgagctggaaaagacctgagtgacaaacagacgcttaccggatcaacgacaaatgccaccactgcgcaacacatcagacggtgtcggctacttccgaggagtgagaggggtcgatgcaactgtctcagtgctgagctctgggccaatttacaga  
tgacaatatctccgctccagagtgaaacgggcagcttccgctaagctaagagacgacattgaaaacctctgtgtgtgactgccaatgagatgtggaatcagttcaacaaagtgaactgtcttcccaatcgattgctgagactgcagatgctaagaataagattcagacgacttagcaa  
agaccctgcagagattttccagactgaaatgacatagaatccatcaagaaggccatcaaggacaagactgcttctgaagtggtcagaccagactggatgagcgcacaagacggccgaacattgagttgtgcccagacatggctcagctacgcttgttaacgaggtacacgag  
gttgacgacacatccagaccctgcagcagcgtgagggatcagaggacaccctgcagtcgctgtgctccatcaaaagccaccctcagatgacctggctgtcaaaagccaattcctgtacatcgaccagaaaaatgatgagcatgcccagagctaccccaacacccctccggc  
tggtcggctctgctag