



Reference: PTXBC031870

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAGK; DAGK1; DGK-alpha; diacylglycerol kinase alpha; 80 kDa diacylglycerol kinase; DAG kinase alpha; diacylglycerol kinase, alpha 80kDa; diglyceride kinase alpha		

For research use only.

GENE DESCRIPTION	diacylglycerol kinase, alpha 80kDa
GENE ID	1606
ORIGIN SPECIES	Human
ACCESSIONS	BC031870
VECTOR	pDONR223
NCBI SYMBOL	DGKA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr