

CSMD2-CUB and Sushi multiple domains 2 Gene

Reference: PTXBC031871

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ1007G16.1; dJ1007G16.2; dJ947L8.1; CUB and sushi domain-containing protein 2; CUB and Sushi (SCR repeat) domain; CUB and sushi multiple domains protein 2; CUB and Sushi multiple domains 2		

For research use only.

GENE DESCRIPTION	CUB and Sushi multiple domains 2
GENE ID	114784
ORIGIN SPECIES	Human
ACCESSIONS	BC031871
VECTOR	pDONR223
NCBI SYMBOL	CSMD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgatgccaccttcgaggccctcgggcatcatcacctcccccaatttccccattcagatgacacaatgcacactgtgtgtggatcatcacagcactcaacccctccaaggtgggaggagaccagggaacatacatgcaacaaggagggtctaagtgtgaaactcaaagtct
aagactgaggaaatagaagttgtcttctgggaaacagcaaggacaatggatgctgtgtgcctaacctgttcatggaaactggcctacaggtgatcaagctcgctttgaggagtttgattggagaggggctatgacacccctgacggctggtgatggtgaggtgaggtgggaccagaagaca
gttctctacatcctgacaggtacatcggtcccgatctcatgtcagcaccaatcaaatgtggctctcttccagactgagtcagtggtgaggttccctgggattcaagcttctatgaagagatcagcagggcaggtgacgtgacctggcatacctgcatalggccggagggaaggct
cccgggttcaccacggtgacacactcaagtttgagtgccagcccgctttgagctggtgggacagaaggcaatcacatgccaaaagaataaccaatgtcggctagaagccaggctgcgtgttaagtga