

FAM71E2-family with sequence similarity 71, member E2 Gene

Reference: PTXBC031875

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM71E2; C19orf16; family with sequence similarity 71 member E2		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 71, member E2
GENE ID	284418
ORIGIN SPECIES	Human
ACCESSIONS	BC031875
VECTOR	pDONR223
NCBI SYMBOL	FAM71E2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatctgcttgaacaggaggtgacctgagccgtccaggcaccgccgaagtgggtccctgtctgggggagctgcaaaagaccctccagaaggcgagtagctgtccctccgtccgtgccatgttcgagagtaactttgtcaggtgacctatcaagggggccagtggtctgtg
aatcacagaaccaaccggctggccatggcgctggccgtccctgccaggcctgtgtgtgacctgacatctgtgatcgccagcccgagggacaggagactgctccgctcgtgtgaccaggatgatccccctggacctgtccacctctggtccatgacctctgctggcg
ctgaagctgcgctgtgtctggccgagtaactgacctggacgccctgacaacgaggtgggtctctgttccactgctgggtccgctcatcaactgcttcaggagccggctccacctggacccccaggaccacgcgcacggccccctggatatgacctggccaaag
cgctgctccacctggcacctgcaggaccagccatcagcagacatgcagtcagtggttgcga