

C7orf61-chromosome 7 open reading frame 61 Gene

Reference: PTXBC031966

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf61; chromosome 7 open reading frame 61		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 61
GENE ID	402573
ORIGIN SPECIES	Human
ACCESSIONS	BC031966
VECTOR	pDONR223
NCBI SYMBOL	C7orf61
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttggatcatgaagttctccgatgggttagacgggttgcaaaaggattattctctgggtttttctctggaggcaaaaattaaaccaaccatctcaggacacccctgactccaagaaacactcattgaagaagatggagaagactctccaggtggtgagactttgaggttggtcagctc
ccaaagaggctaaagcccaagttgggtgagtcctcccgagctggcagatccctgcgtgttgcccaagactacagagagacccaggttgagctgggccaacaggcccaatccctactgcagctgccgagagacggccgtcaagctctctccacgctcatggtctctgccctgcagagc
ggctggcagatgtgagctggaagtcacatgagttctgcctcagtcagctcccaagtgaggacgcagtcaccttgaagactccggaggctgagttgctgtgggaggtgtacctggtgctgtggccgttcggaacacctgcgccggctgtaccgagcaggagagggcacagac
ggcaccacgtccgatgccatgctgcccccgaccaacccggctcagtcacctgaaactggatgcccgaagtcctctag