

C14orf65-chromosome 14 open reading frame 65 Gene

Reference: PTXBC031990

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 65
GENE ID	317762
ORIGIN SPECIES	Human
ACCESSIONS	BC031990
VECTOR	pDONR223
NCBI SYMBOL	C14orf65
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggctcagtgggggccctccacagacagagccgtggccctcaaggtggacagtggagccgctgctcctcagctgcctctgtctcatcacaggccagaagtggaaaggggttcaaatcaagcccccgggtctctgtgtgcgttctacgcagcaaggcttcaagtttcaaaact
gggttggaaactcaaccaggtctcctcaaggaacaaaggcagggcacaggtggggcagggctggcttggcaggcagctgtgcggatccaccagagtgccctgggagggcagtggggttggctgtggctgtgtgagctgcagccaggcctctcccgcgtcccagggtctggaagg
agggcctctgcacatttataagccacagactgcgaaggcctggagccccaccttccaccataacaaagaaaggatagtttcagaactgttttatttggataactgactga