

SPINK6-serine peptidase inhibitor, Kazal type 6 Gene

Reference: PTXBC032003

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BUSI2; UNQ844; serine protease inhibitor Kazal-type 6; acrosin inhibitor; kallikrein inhibitor; protease inhibitor H; serine peptidase inhibitor, Kazal type 6		

For research use only.

GENE DESCRIPTION	serine peptidase inhibitor, Kazal type 6
GENE ID	404203
ORIGIN SPECIES	Human
ACCESSIONS	BC032003
VECTOR	pDONR223
NCBI SYMBOL	SPINK6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactgcaggcatgtttctgctcctctctctgctcttttctgttttaacaggtgtcttcaagtcaggagacaggtgactgtggtgagttccaggacaccaaggtctactgcactcggaaatctaaccacactgtgctctgatggccagacatattgcaataaatgtcctctgtaaggccatagtgaaaagtgtggaagattagcctaaagcatcctgaaaatgctga