



Reference: PTXBC032009

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	glycerophosphoinositol inositolphosphodiesterase GDPD2; GDE3; OBDPF; glycerophosphodiester phosphodiesterase 3; glycerophosphodiester phosphodiesterase domain-containing protein 2; osteoblast differentiation promoting factor; glycerophosphodiester phospho		

GENE DESCRIPTION	glycerophosphodiester phosphodiesterase domain containing 2
GENE ID	54857
ORIGIN SPECIES	Human
ACCESSIONS	BC032009
VECTOR	pDONR223
NCBI SYMBOL	GDPD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]